

TABLE 3. Fast Performance Liquid Chromatography (FPLC) Columns

Column	MW range	Type	Bed volume (ml)	pH range	Capacity (mg)	Ordering information ^a
Mono Q HR 5/5	up to 10^7	anion exchange	1	3–11	25	17-0546-01
Superdex 200 HiLoad 16/60	1×10^4 – 6×10^5	gel filtration	120	3–12	60	17-1069-01
Superdex 200 HR 10/30	1×10^4 – 6×10^5	gel filtration	24	3–12	5–12	17-1088-01
Superose 6 HR 10/30	5×10^3 – 5×10^6	gel filtration	24	3–12	5–12	17-0537-01

^aAmersham Biosciences, www.amershambiosciences.com

Protein–Protein Interactions: A Molecular Cloning Manual, 2nd Ed., © 2005 by Cold Spring Harbor Laboratory Press, Chapter 3, Table 3.